

## Near infrared grain analyzer

**IR2300**

### Features

- Light path design for near-infrared transmission.
- Powerful and professional infrared modeling software.
- The measurement is green and environmentally friendly, no sample preparation is required, and the measurement speed is fast.
- High-definition touch screen operation, intuitive and friendly user experience.
- Long life light source (more than 5000 hours), easy maintenance.
- A variety of data transmission interfaces to facilitate data sharing with other devices.

### Description

IR2300 is a special instrument developed by optosky using the latest near-infrared technology for grain analysis. It delivers superior analytical accuracy and speed, ease of use and robustness. It can conduct non-destructive and rapid detection of moisture, protein, fat, test weight, starch and other parameters in grains.

The instrument adopts touch screen operation, and the software interface is intuitive and concise. Users only need to pour in a grain sample to analyze.



## 1. Principle

Based on the specific spectral characteristics of each substance, the grain sample is irradiated with light of a specific wavelength, and the protein, moisture, and fat in the grain are calculated and analyzed by detecting the light absorption, reflection, or transmission characteristics of the sample, combined with the built-in correction model. , starch and other nutritional components and quality indicators, thereby achieving non-destructive, rapid and accurate detection of grains.

## 2. Application fields

- Purchase of wheat and barley
- Wheat and barley processing
- Wheat and barley sampling testing
- Corn and soybean processing
- Corn and soybean production
- Other cereals, etc.

## 3. Parameter

| Analysis parameters   |   |
|-----------------------|---|
| Samples               | Wheat, barley, soybeans, corn, rice, etc.                     |
| Analysis Time         | About 60 seconds  |
| Parameters            | Moisture, protein, fat, test weight (optional), starch, fiber |
| Principle             | Continuous scanning spectrum, transmission method detection   |
| Sample size           | Whole sample 50-800ml   |
| Number of subsamples  | Customizable  |
| Optical parameters    |   |
| Detector              | 256 InGaAs detector (ATP2000P)                                |
| Wavelength range      | 570-1100nm  |
| Wavelength accuracy   | < 0.05nm  |
| Spectral bandwidth    | About 7nm   |
| Spectral resolution   | 0.1nm~10nm optional   |
| Noise                 | < 20uA  |
| Light source lifespan | Not less than 5000 hours                                      |
| General parameters    |   |
| Interface             | USB, network port   |
| Power supply          | 110V ~ 240V, 50/60Hz  |
| Dimension             | 400*400*350mm   |
| Weight                | 15kg  |