

Portable Raman Spectrometer

ATR3000FD

For Food Safety

Features:

- Over 140 models detected, including drugs, dietary supplement, cosmetics, chemicals;
- Detect limits reach up to ppb level;
- Accurate auto-identification algorithm, no manual readout required;
- High performance, high stable SERS enhanced reagent, life span over 1 year;
- Sample state : Solid, powder, thin film, liquid;
- Touch screen: 12inch;
- Stable Windows 10 operation system;
- USB 2.0, Bluetooth, WIFI, 4G modules (optional);
- Touch screen and keyboard input;
- Friendly interface;
- Battery life> 5hrs;
- IP67;

Description:

Raman spectroscopy is a fast, non-detructive, safety detection technology. It possessed fast, accurate, good repeatabilities. And their samples can be easily pretreated with a compact size multiple functions pretreatment instrument.

ATR3000FD combines SERS with Raman Spectrometer technology, It provides fast identification, in-field detect food pesticide residues, veterinary drugs residues, non-edible chemical substance trace, abuse of food additives, illegal additives drugs for health care.

Models	Range (cm ⁻¹)	Resolution (cm ⁻¹)
ATR3000FD-27	250-2700	4-6
ATR3000FD-35	200-3500	6-8
ATR3000FD-43	200-4300	8-10

Remark:

- Reference to ASTM E2529-06 Method Test ;
- Special customization can improve resolution, but less sensitivity.



Product data information is current as of publication data. Products conform to specifications per the terms of Optosky Standard warranty.

1

Copyright © Optosky Technologies, 2015 Floor 22th, Creative Bld., 1300 Jimei Ave, Jimei, Xiamen, 361005, China Tel: +86-592-6102588



1.1 Specification

ATR3000FD System				
Interface	USB 2.0/LAN			
Touch Screen	12" Touch Screen			
Evidence Picture	3-mega Camera			
Operating System	Windows 10			
Battery life	>5 h			
Integration Time	4ms - 120s			
Power Voltage	DC 19V(+/-5%)			
Operating Temp	-25~60 °C			
Operating Humidity	< 95%			
Storage Condition	-35~75 °C, <95%			
Size (L*W*H)	$35 \times 24.5 \times 15.2 \text{ cm}^3$			
Weight	9 Kg			
Reliability				
Spectral Reliability	$\sigma/\mu < 0.5\%$ (COT 8 hours)			
Temp Reliability	Shift $\leq 1 \text{ cm}^{-1}$ (10-40 °C)			
Spectral Intensity Change	<±5%			
(in 5 ~ 40 °C)				
Optical Parameters				
Spectral Range	250-2700 cm-1 (Customized)			
	250-2700 cm-1 (Customized) 4-6 cm-1 (Customized)			
Spectral Range				
Spectral Range Resolution	4-6 cm-1 (Customized)			
Spectral Range Resolution SNR	4-6 cm-1 (Customized) >570: 1			
Spectral Range Resolution SNR Optical Path	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T			
Spectral Range Resolution SNR Optical Path Focus Distance	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T			
Spectral Range Resolution SNR Optical Path Focus Distance Detector	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels Dynamic Range	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048 1300: 1			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels Dynamic Range Sensitivity	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048 1300: 1			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels Dynamic Range Sensitivity Laser	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048 1300: 1 130 Photon @ 400 nm; 60 Photon @ 600 nm			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels Dynamic Range Sensitivity Laser Central Wavelength	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048 1300: 1 130 Photon @ 400 nm; 60 Photon @ 600 nm 785nm (+/-1nm)			
Spectral Range Resolution SNR Optical Path Focus Distance Detector Model Detector Range Effective Pixels Dynamic Range Sensitivity Laser Central Wavelength	4-6 cm-1 (Customized) >570: 1 f/4, Crossed C-T 98 mm for incidence and output Cooled CCD 200-1100 nm 2048 1300: 1 130 Photon @ 400 nm; 60 Photon @ 600 nm 785nm (+/-1nm) 0.08 nm			

2

Product data information is current as of publication data. Products conform to specifications per the terms of Optosky Standard warranty.



Working Distance	6 mm
Rayleigh Scatter Resistant	OD>8
NA	0.3
Aperture	7mm

1.2 Test Items

1.Application:

- (1) Raman technology is suitable for a wide range of unknown chemical substances;
- (2) The database can be used to measure a wide range of chemical substances including aqueous solution. Sample

shape: constant solid, powder, film, liquid;

(3) System weight ≤ 10 kg, size ≤ 20 cm * 30 cm * 40 cm.

2.Technical Parameter

Portable

Raman

Spectrometer

- (1) Wavelength range: $250-2700 \text{ cm}^{-1}$;
- (2) Resolution: 6cm-1 (Half peak width);
- (3) 2048 pixel CCD detector with refrigeration;

(4) The fully integrated Raman food safety analyzer system based on OEM adopts planar reflection grating spectrometer and Raman laser, and adopts high coupling efficiency and high sensitivity Raman signal coupling and collection;

(5) Weight: ≤ 10 kg,; operation interface: 12 inch industrial true color display screen, multi touch capacitive screen, full Chinese operating system, software can automatically analyze. Automatic analysis of Raman spectrum, intelligent identification of mixture, automatic analysis of mixture composition, without the help of external software and computer processing;

(6) Laser: 785 nm narrow linewidth Raman laser, linewidth better than 0.08 nm, \leq 500 mW,. The laser power is adjustable;

(7) Battery: equipped with 10000 mA lithium battery and lithium battery charging circuit, the battery service time is more than 4 hours;

- (8) Scanning mode: non contact scanning and vial sampling mode;
- (9) Working temperature: 10 °C 40 °C. Working humidity: humidity \leq 90%.
- (10) The circuit board requires Bluetooth module, wireless network card and Unicom 4G communication card;
- (11) Integrated USB interface 2, integrated buzzer;
- (10) Spectrum library: Raman enhanced spectrum of illegal additives in drugs, food and health care products.

3



Users can label more than 50 chemicals in the database. If the detected substance contains the marked chemical substance, the system will display the information containing the prohibited substance. It has the function of adding database items manually;

(11) Data export: report can be exported through SD card;

(12) 24-hour spectrum technical support to help users solve technical problems;

(13) The detection limit reached ppb level, and the detection items included drugs, health care drugs, cosmetics and chemicals;

(14) Sample shape: constant solid, powder, film, liquid;

(15) Detection of illegal additives, agricultural and veterinary hormone residues and adulterated harmful substances in micro trace food;

(16) Detection of illegal additives in micro trace health food;

- (17) Ppb level detection can be achieved by nano enhancement reagent;
- (18) The stability of nano enhanced reagent was more than 12 months.

3. \bigstar Test items (at least, but not limited to the following items, the instrument can expand the test items infinitely):

NO	Item	Common food (scope of application)	Detection time / min	Limit mg/kg
Abuse of food additives				
1	Erythrosine	Drinks, etc	5	0.02
2	Saccharin sodium	Liquor, etc.	5	0.1
3	Aspartame	Drinks, etc	5	0.5
		Non edible chemicals		
4	Urotropine	Bean products such as bean curd sticks	15	0.01
5	Sudan I	Chili oil, chili powder, etc	15	0.2
6	Rhodamine B	Chili oil, chili powder, etc	15	0.2
7	Auramine O	Bean products, yellow croaker, etc	15	0.1
8	Morphine	Seasoning, etc	15	0.1
9	Basic orange 21	Bean products such as bean curd sticks	15	0.1
10	Basic orange 22	Bean products such as bean curd sticks	15	0.1
11	Melamine	Milk and dairy products	10	0.2
12	Sodium thiocyanate	Liquid milk	10	0.05
13	Bengal Rose	Drinks	10	0.1
14	Fuchsia S	Drinks	10	1
15	Sodium fluorescein	Drinks	10	0.1

4



16	Fluorescent Pink B	Drinks	10	1		
17	Patent blue	Drinks	10	5		
Pesticide residue						
18	Malathion	Vegetables, fruits	20	0.1		
19	Fenitrothion	Vegetables, fruits	20	0.5		
20	Parathion methyl	Vegetables, fruits	20	0.1		
21	Parathion	Vegetables, fruits	20	0.2		
22	Trisulfite	Vegetables, fruits	20	0.1		
23	Phorate	Vegetables, fruits	20	0.2		
24	Phoxim	Vegetables, fruits	20	0.2		
25	Thiophos	Vegetables, fruits	20	0.25		
26	Methamidophos	Vegetables, fruits	20	0.1		
27	Acephate	Vegetables, fruits	20	0.1		
28	Imiphos	Vegetables, fruits	20	0.5		
29	Fenthion	Vegetables, fruits	20	0.05		
30	Chlorpyrifos	Vegetables, fruits	20	0.5		
31	Chlorpyrifos methyl	Vegetables, fruits	20	0.5		
32	Triazophos	Vegetables, fruits	20	0.5		
33	Isocarbophos	Vegetables, fruits	20	0.5		
34	Dicofol	Vegetables, fruits	20	1		
35	Difenoconazole	Vegetables, fruits	20	0.1		
36	Pyridaben	Vegetables, fruits	20	0.1		
37	Carbaryl	Vegetables, fruits	20	1		
38	Paraquat	Vegetables, fruits	20	0.01		
39	Triadimefon	Vegetables, fruits	20	0.1		
Veterinary drug residues						
40	Malachite green	Aquatic products, fish	15	0.01		
41	Hidden malachite green	Aquatic products, fish	15	0.01		
42	Enrofloxacin	Aquatic products, fish	15	0.02		
43	Sodium pentachlorophenol	Aquatic products, fish	10	5		
44	Crystalloid violet	Aquatic products, fish	10	0.01		
	-	Illegally added in health products				
45	Sibutramine	Slimming health care drugs	10	0.5		
46	Ephedrine	Slimming health care drugs	10	0.5		
47	Sildenafil	Anti fatigue health care drugs	10	0.1		
48	Tadalafil	Anti fatigue health care drugs	10	0.1		
49	Edenafil citrate	Anti fatigue health care drugs	10	0.5		
 50	Reddenafil	Anti fatigue health care drugs	10	1		

5



51	Vardenafil	Anti fatigue health care drugs	10	0.05
52	Pseudo Vardenafil	Anti fatigue health care drugs	10	1
53	Triazolam	Tranquilizing and sleeping care drugs	10	0.1
54	Diazepam	Tranquilizing and sleeping care drugs	10	0.1
55	Oxazepam	Tranquilizing and sleeping care drugs	10	0.1
56	Clozapine	Tranquilizing and sleeping care drugs	10	0.1
57	Clonazepam	Tranquilizing and sleeping care drugs	10	0.1
58	Alprazolam	Tranquilizing and sleeping care drugs	10	0.1
59	Captopril	Antihypertensive health care drugs	10	0.5
60	Clonidine hydrochloride	Antihypertensive health care drugs	10	0.5
61	Prazosin hydrochloride	Antihypertensive health care drugs	10	0.5
62	Nifedipine	Antihypertensive health care drugs	10	0.5
63	Atenolol	Antihypertensive health care drugs	10	0.1
64	Glipizide	Hypoglycemic health care drugs	10	1
65	Gliclazide	Hypoglycemic health care drugs	10	0.5
66	Glibenclamide	Hypoglycemic health care drugs	10	0.5
67	Repaglinide	Hypoglycemic health care drugs	10	1
68	Metformin hydrochloride	Hypoglycemic health care drugs	10	0.1
69	Pioglitazone hydrochloride	Hypoglycemic health care drugs	10	0.5
70	Phenformin hydrochloride	Hypoglycemic health care drugs	10	0.1
71	Rosiglitazone maleate	Hypoglycemic health care drugs	10	0.1

1.3Other Accessories

For the first order, 200 times of SERS enhancement reagents were purchased as gifts;





1.4Successful Application Cases

Time	Customer	Item	Application
August 2017	Weifang Public Security	Raman integrated rapid screening	Rapid screening of food safety
	Bureau	equipment for food safety	
September 2017	Haizhou District,	Raman integrated rapid screening	Rapid screening of food safety
	Lianyungang	system for food safety	
December 2017	Dezhou City, Shandong	Raman integrated food safety	Rapid screening of food safety
	Province	screening system	
December 2017	Sichuan food and Drug	Food safety emergency rapid	Rapid screening of food safety
	Administration	inspection equipment	

1.5Optional Accessories

***Optional:** The ATFD030 multi-functional pretreatment equipment, which is specially used for **surface enhanced Raman detection, integrates the functions of ultrasonic extraction, centrifugation and volatile concentration**. It is portable, automatic and easy to operate. It can be used for the separation and enrichment of complex system in the process of surface enhanced Raman detection of actual samples.



Fig. 1 Physical diagram of ATFD030 (vehicle mounted) portable multifunctional pretreatment instrument