

Confocal Raman Microscope

ATRH8800

Features

- Microscopic hyperspectral imaging, confocal micro-Raman spectroscopy, microscopic imaging;
- Micro-area positioning device, laser spot alignment and adjustment
- High power optical positioning system
- Automatic Raman imaging experiment, automatic focusing, automatic scanning;
- Wide-band range hyperspectral imaging (390-1000nm);
- Automated hyperspectral imaging experiment;
- SNR Ratio >6000:1
- Unique software controls switching light path
- Quickly locate focus position
- High-quality objective lens 5 million cameras
- Excitation wavelength: 532, 633, 785, 830, 1064nm optional
- USB2.0

Application

- Nanoparticles and new materials
- Universities and research institutes
- Biology
- Forensic Medicine Identification
- material science
- Medical Immunoassay
- Agriculture and food identification
- water pollution analysis
- Gem and inorganic mineral identification

Description

ATRH8800 adds a microscopic hyperspectral imager to the confocal Raman microscope. The instrument is an advanced device that integrates a hyperspectral imaging microscope, an optical microscope and a high-resolution confocal Raman microscope. Hyperspectral imaging microscopy and confocal micro-Raman spectrometer can be used to characterize and analyze the surface morphology, reflection hyperspectral imaging and Raman spectrum performance of nanomaterials respectively, thereby providing more comprehensive information on the sample and providing sharp Microscopic image. Users can improve work efficiency, spend more time on data collection and analysis, and truly realize in-situ detection and analysis of samples.

ATRH8800 is equipped with an objective lens specially designed for the Raman system, which makes the laser spot close to the diffraction limit, and then displays the focus information accurately and intuitively on the computer through a 5-megapixel camera. It overcomes the problem in ordinary Raman systems that the focal plane for collecting Raman signals is slightly higher or slightly lower than the actual optimal focal plane, thereby improving the quality of the Raman spectrum.

At the same time, ATRH8800 uses high-performance Raman specially optimized for micro-Raman systems. It is industry-leading in terms of sensitivity, signal-to-noise ratio, stability, etc., providing a strong guarantee for Raman research.



1. Parameter

ATRH8800 Parameter	
Excitation Lasers	Select up to four from 266、325、532、638、785、1064nm
Laser Power	266nm: 30 mW 325nm: 30mW 532nm: 100mW 633nm: 80mW 638nm: 80mW 785nm: 350mW 1064nm: 500mW
Spectroscopic Spectrum	Asymmetric C-T optical path
Focal Length	350mm、510mm、810mm optional
Built-in Gratings	Standard configuration: 3 pieces; 300 lines, 600 lines, 1200 lines, 1800 lines, 2400 lines optional
Detector	1) Deep cooling area array CCD: 2000X256 pixels
Objective Lens	2) Deep cooling high sensitivity EMCCD: 1600X200 pixels
Microscopic illumination	3) Deep cooling area array InGaAs CCD: 512X1 pixels
Lighting Method	Up to 2 detectors can be integrated, choose one from detector 1# or detector 2#;
Camera System	Standard configuration: 4X, 10X, 20X, 50X;
Focus Method	Optional configuration: 100X
Laser Spot Diameter	High brightness and long life white LED
Laser Stability	epi-illumination
Communication Mode	5 million pixel industrial camera
Hyperspectral imaging parameters	
Imaging Method	Push broom
Spectral Range	390- 1000nm
Spatial Channels	1200
Wavelength Channels	1920
Spectral Resolution	1.3 nm
Imaging Speed	<60s
X, Y axis two-dimensional platform	
Movement Method	Fully electric, manual optional
Imaging Range	50 X 50 mm, 100 X 100 mm optional
Mobile Resolution	0.1 μ m
Positioning Accuracy	1 μ m
Scan Interval	Software settings, minimum 1 μ m

Scan Speed	20mm/s
Nanodisplacement Stage (Optional)	Minimum displacement resolution 2nm, displacement accuracy 10nm
Z axis (auto focus)	
Focus Accuracy	$\leq \pm 0.2 \mu\text{m}$
Maximum Stroke	20 mm
Focus Speed	Less than 10 s
Nanodisplacement Stage (Optional)	Minimum displacement resolution 2nm, displacement accuracy 10nm
Physical parameters	
Size /mm	ATRH8800-FL350: 905(L)×58.3(W)×643(H) ATRH8800-FL510: 1009(L)×58.3(W)×643(H) ATRH8800-FL810: 1520(L)×68.3(W)×643(H)
Weight	ATRH8800-FL350: 59 Kg ATRH8800-FL510: 63 Kg ATRH8800-FL810: 78 Kg
Working environment	
Voltage	100~240VAC
Peak power	<500W
Other motivations	No need
emission	none
Platform requirements	Air-floating vibration isolation optical platform
Working temperature and humidity	Constant temperature ($25 \pm 2^\circ\text{C}$), constant humidity ($50 \pm 10\%$)
cleanliness	Level 10,000 and above



Figure 2 ATRH8800 microscopic Raman functional structure indication diagram

Product data information is current as of publication data. Products conform to specifications per the terms of Optosky Standard warranty.




Figure 3 Software interface of ATRH8800






Figure 4 Software interface of ATRH8800






Figure 5 test results indicate: the instrument resolution reaches 0.847 cm-1, test specification: tested according to the national standard "General Specifications for Raman Spectrometers" method, test equipment: ATRH8800-FL510, test light source: mercury argon lamp, collection spectrum line: 546.08 nm




Figure 6 Spectrum of ATRH8800-FL510 test

2. Selection Guide

Table 1 ATRH8800 Selection Guide

Model	Focal Length	excitation wavelength/nm*3	Laser power/mW	Maximum wave number range*1*2	Optimal resolution/cm-1*4
ATRH8800-FL350	350mm	325	30	50~ 10000	2.6
		532	100	5 ~ 10000	1.4
		638	80	5 ~ 10000	1.4
		785	350	5 ~ 10000	2.1

Product data information is current as of publication data. Products conform to specifications per the terms of Optosky Standard warranty.

		1064	500	50 ~ 10000	5.2
ATR8800-FL510	510mm	325	30	50~ 10000	1.9
		532	100	5 ~ 10000	0.9
		638	80	5 ~ 10000	0.9
		785	350	5 ~ 10000	1.4
		1064	500	50 ~ 10000	3.6
ATR8800-FL810	810mm	325	30	50~ 10000	1.1
		532	100	5 ~ 10000	0.45
		638	80	5 ~ 10000	0.45
		785	350	5 ~ 10000	0.86
		1064	500	5 ~ 10000	2.3
ATR8800LT: deep cooling to -30°C, ultra-long integration time (up to 1.3h) ATR8800EM: Deeply cooled area array EMCCD detector ATR8800BS: basic type ATR8800AF: Autofocus type ATR8800MP: Scanning imaging Mapping, autofocus					

Note:

*1: Maximum wave number range, related to the selection of excitation wavelength;

*2: The standard configuration is 150cm-1, and the minimum wave number range can be customized to 5cm-1, 50cm-1;

*3: Other excitation wavelengths can be customized;

*4: The optimal resolution is related to the slit width of the spectrometer; if the slit width is further reduced, the resolution can be further improved;

*5: The parameters in the table only represent the parameters of the company's standard products; Aopu Tiancheng instruments are all independently developed and produced products, and the corresponding parameters can be customized;

Ordering Guide:

Naming example:

ATR8800AF-LT-FL350-532+638: Autofocus, long integration time, focal length of 350mm, dual excitation wavelengths: 532nm and 633nm respectively

ATR8800MP-EM-FL810-532+638+1064: scanning imaging, EMCCD detector, focal length of 810mm, excitation wavelength of three wavelengths: 532nm, 633nm and 1064nm respectively