

ATP5030/4

Cooled, Ultra High Resolution

Features

- Detector cooling temperature: -5 °C;
- M type optical path, higher resolution;
- Spectral range: 200-1100 nm;
- Spectral resolution: 0.05 ~ 2 nm;

Micro Spectrometer

- Detector: 2048 or 4096 pixels;
- Integration time: 0.1ms ~ 256s;
- Power supply: DC 5V power supply;
- ADC bit depth: 16 bits;
- Data output interface: USB Type-C;
- 20-pin expansion interface.

Applications

- LIBS, Plasma luminescence detection;
- Raman spectroscopy detection;
- Wavelength monitoring, laser, LED and other luminous bodies;
- Water quality analyzer;
- UV flue gas analyzer;
- LED sorting machine, color detection;
- Spectral analysis, radiation spectroscopic analysis, spectrophotometric analysis;
- Reflection and transmission spectrum detection.

Description

ATP5030/4 is a cooled, ultra-high resolution micro spectrometer newly developed by Optosky, ATP5030/4 is based on the M-type C - T optical path structure independently developed by Optosky, with extraordinary ultra-high resolution; at the same Optical path structure: non-crossing M-shape C-T time, it uses a cooled 2048 or 4096 pixel linear array detector to achieve the ultimate ultra-high resolution, the highest resolution can reach 0.05nm, It is a miniature spectrometer suitable for various high-resolution applications, with high reliability, ultra-high speed, low cost, and high cost performance. It can be adapted to various environmental uses such as online testing.


> ATP5030/4 adopts TEC electric cooling, the cooling temperature is -5°C. Greatly reduce the dark current and noise of the spectrometer, and improve the dynamic range and signal-to-noise ratio of the spectrometer.

> ATP5030 can receive SMA905 optical fiber input light or free space light, and output the measured spectrum data through USB2.0 or UART port.

> ATP5030 only needs a 5V DC power supply, which is very easy to integrate and use.

	Detector pixels	
ATP5030	2048	
ATP5034	4096	
ATP5030P	2048	Cooled,
AIPSUSUP	20 4 8	Back-thinned CCD





Product data information is current as of publication data. Products conform to specifications per the terms of Optosky Standard warranty.




1.Parameter


Detector	
Туре	Linear array detector
Detectable range	200-1100 nm
Effective pixel	2048 or 4096 pixels
Sensor Cooled	TEC cooled, -5 °C
Pixel dimension	14μm × 200μm
Sensitivity	1300 V/(lx·s)
Dark noise	13 RMS @ 13 °C
Optical Parameter	
Wavelength range	200-1100 nm
Optical resolution	0.07-3 nm
Signal-to-noise	>600:1
Dynamic range	8.5 x 10 ⁷ (system); 2000:1 for a single acquisition
Stray light	<0.05% at 600 nm; <0.09% at 435 nm
Optical Configuration	
Optical Design	M-type C-T
Focal Distance	75mm
Incidence slit	50 μm (10, 25, 100, 200 μm are optional)
Incident Interface	SMA905 connector
Electrical Parameter	
Integration time	0.1 ms - 256 second
Interfaces	USB Type-C
A/D conversion resolution	16 bit
Supply voltage	DC4.5 to 5.5 V (type @5V)
Operating current	1.5 A@Typ. 3A Max
Storage temperature	-30°C to +70°C
Operating temperature	-25 ~ 50 °C
Working humidity	< 90%RH
Physics Parameter	
Dimension	135 × 99 × 46 mm
weight	0.5 kg
Sealing	Anti-sweat




2 Mechanical Diagrams
























3 Electrical Pin-out

Table 1 Electrical Characteristics

Parameter	Min	Тур	Max	Unit
Power Supply				
Operating voltage range	4.5	5	5.5	V
Operating current		170		mA
Logic Inputs(3.3V LVTTL,				
Five-volt tolerant)				
High level input voltage	1.7		3.6	V
Low level input voltage	-0.3		1.0	V
Logic Output(3.3V LVTTL)				
High level output voltage	2.4			V
Low level output voltage			0.4	V

The module is equipped with a 30-pin male angled box header(2x15, 2.00 mm pitch) and Type-C interface.

Table 2 Electrical Pin-Out

Pin#	Description	I/O	Function Description		
1	VCC	1	Power Supply, 5V \pm 0.5,		
2	GND	1	Ground		
3	UART_TX	Output	UART Transmit signal		
4	UART_RX	Input	UART Receive signal		
5	Lamp_En	Output	LVTTL output the lamp enable signal.		
6	Continuous_ strobe	Output	LVTTL output the continues strobe signal.		
7	Ext_trigger_i	Input	LVTTL input the trigger signal.		
8	Single_strob	Output	LVTTL output the single strobe signal.		
9	SPI_SCK	Output	The SPI Clock signal for communications to other SPI peripherals		
10	SPI_MOSI	Output	The SPI Master Out Slave In (MOSI) signal for communications to other SPI peripherals		
11	SPI_MISO	Input	The SPI Master In Slave Out (MISO) signal for communications to other SPI peripherals		
12	SPI_CS	Output	The SPI Chip/Device Select signal for communications to other SPI peripherals		
13	GPIO0	Input /Output	General Purpose Software Programmable Digital Inputs/Outputs, LVTTL Logic.		



1.1	CDIO4	Input	General	Purpose	Software	Programmable	Digital
14	GPIO1	/Output	Inputs/Outputs, LVTTL Logic.				
15	5 GPIO2	Input	General	Purpose	Software	Programmable	Digital
13	GFIOZ	/Output	Inputs/Outputs, LVTTL Logic.				
16	GPIO3	Input	General	Purpose	Software	Programmable	Digital
10	GPIOS	/Output	Inputs/Ou	Inputs/Outputs, LVTTL Logic.			
17	GPIO4	Input	General	Purpose	Software	Programmable	Digital
17	GPIO4	/Output	Inputs/Outputs, LVTTL Logic.				
10	CDIOE	Input	General	Purpose	Software	Programmable	Digital
18	GPIO5	/Output	Inputs/Outputs, LVTTL Logic.				
19	GPIO6	Input	General	Purpose	Software	Programmable	Digital
19	GPIO0	/Output	Inputs/Outputs, LVTTL Logic.				
20	00 0007	Input	General	Purpose	Software	Programmable	Digital
20	20 GPIO7		Inputs/Outputs, LVTTL Logic.				

4 Order Guide

Order number Rules:

Model	Spectral region		Slit width	
ATP5030	Short wavelength	Long wavelength	Slit width	

For example:

What to buy ATP5030, spectral region: 200-1000nm, slit width is 50 um, then the order no is:

ATP5030-200-1000-050

Order No	Spectral region	Slit
ATP5030-200-400-###	200~400	10 μm
ATP5030-200-850-###	200~850	25 μm
ATP5030-200-1000-###	200~1000	50 μm
ATP5030-340-850-###	340~850	100 μm
ATP5030-600-1100-###	600~1100	200 μm
ATP5030-###-###	Other	Other:µm



5. Company Profile

Optosky company is a first-class spectroscopy solution provider, with the headquarter locates in the 7th floor of the research institute of the Chinese Academic of Science at an area of 2500 square meter in Xiamen city where successfully held the international 9th BRICK summit in 2017. The subsidiary company locates in Wuhu city with an area of 2035 square meters.

The company founder Dr.Hongfei,Liu graduated Docter degree from the Chinese Academic of Science and postdoctoral degree from Xiamen University, by integrating both of top Universities' spectroscopy technology background into Optosky company aiming at developing the leading spectroscopy equipment in the world.

The company bases on unique technologies of Optomechatronics, Spectroscopy Analysis, Process Weak Optical and Electrical Signals, Cloud Computing, and have been developed wide products line of the competitive Raman spectroscopy instruments, micro spectrometer, hyperspectral imager, field spectroradiometer, fluorescence spectroscopy, LIBS etc. Driven by advanced technologies and products, Optosky brand has been well-known to customers all over the world.

Optosky company base on technology innovation, market-driven direction, customer first, provides first-class products and services, and one-stop solutions to many fortune 500 companies in many industries. The company received praise from different industry companies, as well as many innovative intellectual properties, software copyright, qualification certification, and winner awards over hundred numbers.

Optosky receives top class A introduced the high-tech company to international Xiamen city, the national high-tech and new innovative technology company award. The founder Dr.Hongfei Liu receives the innovation talent award by the ministry of science and technology.

The company is currently conducting the exclusive project of major industrialization national oceanic administration with a total fund of five million us dollars. The company in charge of drafting national industry standard of VNIR and SWNIR Field Spectroradiometer, and six national standard drafters, including China National Standard Drafter for Hazmat detector based on Raman spectroscopy, China National Standard Drafter for Buoy-type Monitor eco-environment, China National Standard Drafter for water quality monitor in the unmanned boat, China National Standards drafter for online water quality monitor by spectroscopy, China National Standard Drafter for UV-absorbent measure fabrics.

The company has over 70 IPs and over 20 innovative patents.



The company received ISO9001:2015 certification, CE certification, Police Administration Certification, FDA approval compliant, IQOQPQ compliant.



Figure 1 Optosky (Xiamen) Photonics Inc. Company Headquarter




Figure 2 Optosky Company Area






Figure 3 Oversea Market Shares



Figure 4 Optosky Chair and Draft National Standards Lists.

-9-





Figure 5 Qualification

Informationization & Industrilization Fusion Management System




Figure 6 GB/T 23001 Informationization & Industrilization Fusion Management System





Figure 7 Optosky's Co-founder_Dr. Hongfei Liu




Figure 8 Category & Application






Figure 9 Model Name Rule